



New: ArtiSential Large Clip Applier

The World's First and Only 90° Fully-Articulating, Double-Jointed Clip Applier



Developed for Robotic Surgeons, including:

- Urologists performing Robotic Prostatectomy and Robotic Nephrectomy
- Surgical Oncologists and Hepatobiliary Surgeons performing Robotic and Laparoscopic procedures
- General Surgeons requiring articulation during Robotic and Laparoscopic procedures
- Pediatric Surgeons performing minimally invasive procedures

The **ArtiSential™ Clip Applier** offers an unmatched range of motion with 90° articulation, enabling surgeons to achieve precise angles of approach and flexible maneuverability — even in the narrow pelvic cavity. This advanced flexibility enhances the ability to target critical anatomical structures with greater accuracy supporting more precise clip application in both robotic and laparoscopic surgery.

- Uniquely suited to the anatomy of the pelvis and retroperitoneum
- 6-Degrees of Freedom: Double-wristed, double-jointed design allows for precise targeting of anatomical structures, reducing tissue manipulation for proper clip positioning
- Reduced manipulation may decrease bleeding and the risk of tissue damage
- Improved clamp efficiency: May reduce procedure times by facilitating easier, faster, and more precise clip placement
- Improved jaw strength with our large clip handles larger tissue bundles and structures, saving time
- Can be placed through an 8 – 12mm trocar to assist at the bedside during Robotic procedures



ACA02-LV Clip Applier Large-38cm-Lock Adjustable
ACA02S-LV Clip Applier Large-25cm-Lock Adjustable
ACA01-LV Clip Applier Medium Large-38cm-Lock Adjustable
ACA01S-LV Clip Applier Medium Large-25cm-Lock Adjustable

www.LivsMed.us (619) 872-5555 infoUSA@livsmed.com